

# MedaSonics TRIA® II Handheld Dopplers

## For Greater Dependability, Durability and Accuracy

The all-new MedaSonics TRIA II handheld Doppler utilizes improved SSQ (Superior Sound Quality) and SSP (Sensitec Signal Processing™) electronics for outstanding accuracy and reliable readings. Available with and without an integrated digital Accu-Rate™ information display, the TRIA II comes with additional features including:

- Interchangeable probe capability – Enhances system flexibility
- User replaceable cable – Reduces downtime and costly repairs
- Extensive probe options – Including fetal, waterproof fetal, vascular and sterilizable probes
- Enhanced probe sensitivity – Detect fetal heart rate faster and sooner
- Extended battery life – Lasts longer between battery replacements



# MedaSonics TRIA® II Handheld Dopplers

## Tria II Specifications

- Fetal Doppler: 2.1 MHz and 3.2 MHz
- Waterproof Fetal Doppler: 2.1 MHz
- Vascular Doppler: 4, 5 and 8 MHz (not for fetal use)

## Outputs

One 3.5 mm jack to connect to a headset or audio recording device

## Controls

Separate volume control and push ON/OFF button

## Power Source

- Battery voltage:
  - (3) AA Alkaline 1.5 volt (non-rechargeable)
  - (3) AA NMHi 1.2 volt (rechargeable)
- Battery life:
  - 1,250 one minute exams (Alkaline)
  - 1,000 one minute exams (NMHi) rechargeable

## Weight

11.29 oz. (320 gm)

## Dimensions

5.5 x 2.75 x 1.3 in (140 x 70 x 35 mm)

## Service

CooperSurgical, Inc. maintains a full-service facility that can promptly repair products returned to the factory.

## Return/Warranty

- CooperSurgical, Inc. is pleased to present our customers with a 30-day no hassle evaluation return policy in the event you are not satisfied with the product.
- Dopplers are covered by a 5-year warranty on parts, 1-year warranty on labor.

# TRIA® II Ordering Information

Code	Description
<b>MedaSonics TRIA II</b>	
T334	Non-display, Interchangeable 3 MHz probe
T345	Non-display, Interchangeable 2.1 MHz probe
T334-D	Display, Interchangeable 3 MHz probe
T345-D	Display, Interchangeable 2.1 MHz probe
T335	Non-display, Interchangeable 5 MHz probe
T335-D	Display, Interchangeable 5 MHz probe
T336	Non-display, Dedicated 2.1 MHz Waterproof probe
T336-D	Display, Dedicated 2.1 MHz Waterproof probe
T337	Non-display, Interchangeable 8 MHz probe
T337-D	Display, Interchangeable 8 MHz probe
T337-S	Non-display, Interchangeable 8 MHz sterilizable probe
T337-DS	Display, Interchangeable 8 MHz probe

<b>MedaSonics TRIA II Probes</b>	
TSP-2	Fetal Probe 2 MHz
TSP-3	Fetal Probe 3 MHz
TSP-5	Vascular Probe 5 MHz
TSP-8	Vascular Probe 8 MHz
TSP-2WP	Waterproof 2 MHz Probe
TSP-8S	Sterilizable Vascular Probe 8 MHz

<b>MedaSonics TRIA II Accessories</b>	
CRT0007	Wall Mount
CTCBL0001	Interconnect Cable
CT200	Rollstand with Basket
CT230	Rollstand with no Basket
CTGEL0001	Ultrasound Gel, 60 gm
CTGEL0002	Ultrasound Gel, 250 gm



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

800.243.2974 | 203.601.5200 | [coopersurgical.com](http://coopersurgical.com) | 95 Corporate Drive | Trumbull, CT 06611  
Also on [YouTube.com/coopersurgical](https://www.youtube.com/coopersurgical)